



Installing courier and selecting a Lab


1. First go to <https://www.360courier.net/>
2. Select download.




3. The practice will fill out all the fields and select “register and download”.



4. 360courier installer will download and the user will select “run” and next thru the prompts.



5. The end user will get an automatically generated email with the 360courier serial and activation. The codes can both be copied and pasted right into the 360courier field.



Installing courier and selecting a Lab

This is a copy of the activation email.

To: [REDACTED]

Dear 360imaging,

We are delighted for your interest in our valuable application 360Courier.

360Courier application is there to make it easy for you to share your DICOM images, 360dps plans and other attached files in a HIPPA compliant way, with your lab and also with your colleagues, submit, track and manage your implant planning and surgical guide orders.

360Courier is a valuable addition that complements the wonderful planning experience of your 360dps application.

Please download the latest version of courier from <http://www.360courier.net>

After installation, use the following serial and activation codes to start enjoying your experience.

Serial: 10ec46973adf430a8878023194c43 [REDACTED]
Activation: b778c1567bbb4303a310641cb86f4 [REDACTED]

Please be noted that your service contract is valid till: 2100/01/01

In case of lost or non working serials or any inquiries please inform 360Courier support team on: support@360courier.net


also, it is important to make sure that your personal information is correct as recorded in our database.

Your personal information are:


Name: 360imaging support
Email: [REDACTED]
Mobile: 404-236-7700

Steps on how a doctor selects his Lab.

1. Select the name on the account in the top right of courier.




2. Select preferences.




Installing courier and selecting a Lab

3. Select the Lab.




4. Now select your Lab and select OK.



5. Now select Practice information and fill out the entire form. Save changes

Installing courier and selecting a Lab



The image shows a 'User Preferences' dialog box with three tabs: 'Basic Information', 'Practice Information', and 'Preferences'. The 'Practice Information' tab is selected and highlighted with a red border. Inside this tab, there are two main sections, both also highlighted with red borders. The first section contains a text field with 'TEST' and a dropdown menu with 'Orthodontist'. The second section, titled 'Address', contains four input fields: 'Address Line 1', 'Address Line 2', a dropdown menu with 'USA', and a dropdown menu with 'Georgia'. Below these are two more text fields: 'Dallas' and '30132'. At the bottom right of the dialog box are 'Cancel' and 'Save Changes' buttons.

Now your lab is selected and all files uploaded will directly go to your selected Lab.